

MODEL Higi



A free and easy way to start keeping tabs on your health and wellness information.

Standard Tests

- Blood Pressure
- Body composition analysis (body fat percentage)
- Weight
- Body Mass Index
- Heart Rate
- Hydration Index

Standard Configuration

- 10.4" LCD touchscreen Privacy Screen
- 19" HD LCD display screen
- FDA-cleared components
 - blood pressure cuff (adaptive positioning)
 - body composition analyzer (bioimpedance)
- Rugged scale (max. 440 lbs. / 200 kg.)
- Rugged seat, wheelchair accessible (ADA-compliant)
- Secure internet connectivity / LAN & Wifi Enabled
- Connect over 50 activity tracking devices to Higi using your smart phone or log in to your Higi account at higi.com/login & connect through "Apps & Devices"

Offered

- For Sale
- For Lease

Available

- Local Service



WE GUARANTEE THE BEST PRICING & SUPPORT SERVICES

Ask about our **LOWEST PRICE GUARANTEE** program



MODEL Higi



Custom Features

- Your custom logo and branding
- Dynamic messaging, incentives and advertising
- Secure data hosting and software updates
- Installation, training and service
- Online health monitoring tools
- Comprehensive custom reporting

Operating Temperature

41-95 °F (5 to 35 °C)

Operating Humidity

10-80% (non-condensing)

Service & Compliance

- FDA-cleared class II medical devices
- American Disabilities Act (ADA) Compliant
- HIPAA compliant
- RoHS compliant
- UL/IEC 60601-1 1/2 certified
- IEC 8100-1:2002 compliant blood pressure device
- CAN/CSA-C22.2 certified
- Licensing agreement includes: software, updates, parts, service and standard usage reports

Device Options

- Cellular connectivity hardware



Physical Dimensions

- Dimensions: 65"H x 30.6"W x 33"D / 58.8" w/ wheelchair access
- Weight 350lbs. (158.6 kg.)

