

MODEL LC500



A full automatic self-attended health kiosk offering a comprehensive number of biometric measurements and advance features.

STANDARD TESTS / FEATURES

- Blood Pressure
- Body Mass Index (BMI) Calculator
- Weight
- Heart Rate
- Select Pedometers USB downloadable
- Select Glucose Meters Infra-red downloadable
- Bilingual English/Spanish
- 10.4" LCD Touchscreen
- Interactive Health Information Screens
- Color Code4d Results Charts
- Four Information Display Panels
- Banner Display Panel

DEVICE OPTIONS

- Body Composition Analysis
- Pulse –OX Provides real-time evaluation of blood oxygen levels
- HRA/Health Risk Appraisal (not savable)
- Wheelchair Accessibility



* kiosks image configured with options



WE GUARANTEE THE BEST PRICING & SUPPORT SERVICES

Ask about our **LOWEST PRICE GUARANTEE** program



MODEL LC500



OFFERED

- For Sale
- For Lease

INCLUDES

- 1 year Warranty

AVAILABLE

- Local Service
- Extended Warranty

EQUIPMENT

- Class II medical device manufactured in FDA-regulated facility

ACCURACY

- Complies with AAMI (Association for the Advancement of Medical Instrumentation) standards

TEMPERATURE

- Operating: 50 to 102°F (10 to 39°C)
- Storage: 40 to 131°F (4 to 55°C)

HUMIDITY

- Operating: 10 to 95% (noncondensing)
- Storage: 10 to 95% (noncondensing)

DIGITAL DISPLAY

- Systolic & diastolic pressure, pulse rate

RATE OF CUFF INFLATION

- Approximately 8 mmHg / sec

MEASUREMENT RANGE

- Systolic: +60 to +260 mmHg
- Diastolic: + 25 to +240 mmHg
- Pulse: 40 bpm to 199 bpm

MEASUREMENT METHOD

- Oscillometry

OVERPRESSURE LIMIT

- 305 +/- mmHg

TEST TIME

- Typically less than one minute

INPUT POWER REQUIREMENTS

- 117 VAC, 60 Hz
- 2.5 amps
- 3-wire grounded outlet



* kiosks image configured with options

PHYSICAL DIMENSIONS

- Height 68in. (172.7cm.)
- Width 35.5in. (85cm.)
- Depth 25in. (63.5cm.)
- Weight 175lbs. (79.3kg.)

