

# MODEL 90550 / LC300



## Automated blood pressure and heart-rate monitoring in a space saving design

### #1

### WORLDS MOST WIDELY USED HEALTH KIOSK

#### Features

- Blood Pressure and Pulse Rate
- Proven Accuracy
- Fully Automated Upper Arm Pressure Cuff
- Led Display of Readings
- Blood Pressure Guidelines Display
- Oscillometric Methodology
- Four Information Display Panels
- User Record Card Holder

#### Offered

- For Sale
- For Lease

#### Includes

- 1 Year Warranty

#### Available

- Local Service
- Extended Warranty



**BEST  
PRICE**

**GUARANTEE**

## WE GUARANTEE THE BEST PRICING & SUPPORT SERVICES

Ask about our **LOWEST PRICE GUARANTEE** program



# MODEL 90550 / LC300



## Equipment

- Class II medical device manufactured in FDA-regulated facility

## Accuracy

- Complies with AAMI (Association for the Advancement of Medical Instrumentation) Standards

## Temperature

- Operating: 50 to 102° F (10 to 39° C)
- Storage: 40 to 131° F (4 to 55° C)

## Humidity

- Operating: 10 to 95% (noncondensing)
- Storage: 10 to 95% (noncondensing)

## Digital Display

- Systolic & diastolic pressure, pulse rate

## Rate of Cuff Inflation

- Approximately 8 mmHg/sec

## Measurement Range

- Systolic: +60 to +260 mmHg
- Diastolic: +25 to +240 mmHg
- Pulse: 40 bpm to 199 bpm

## Measurement Method

- Oscillometry

## Overpressure Limit

- 305 +/- mmHg

## Test Time

- Typically less than one minute

## Input Power Requirements

- 117 VAC, 60 Hz
- 2.5 amps
- 3-wire grounded outlet



## Physical Dimensions

- Height 64.3 in. (163.2 cm.)
- Width 24.2 in. (61.5 cm.)
- Depth 20.3 in. (52.6 cm.)
- Weight 90 lbs. (40.1 kg.)

