

# Model 90550 /LC300

Automated blood pressure & heart-rate monitoring in a space saving design.



## FEATURES

- Blood pressure and pulse rate
- Fully automated upper arm pressure cuff
- LED display of readings
- Blood pressure guidelines display
- Four information display panels
- User Record Card molder

### OFFERED

- ✓ For Sale
- ✓ For Lease

### INCLUDED

- ✓ 1 year warranty

### AVAILABLE

- + Local service
- + Extended warranty

Ask about our  
**Best Price Guarantee**



(800) 800.1787

AmericanMedicalScreening.com



## Model **90550 /LC300**

### **DIMENSIONS**

- Height 643 in. (163.2 cm.)
- Width 24.2 in. (61.5 cm.)
- Depth 20.3 in. (52.6 cm.)
- Weight 90 lbs. (40.1 kg.)

### **TEMPERATURE**

- Operating: 50 to 102° F (10 to 39° C)
- Storage: 40 to 131° F (4 to 55° C)

### **HUMIDITY**

- Operating: 10 to 95% (noncondensing)
- Storage: 10 to 95% (noncondensing)

### **CUFF INFLATION RATE**

- Approximately 8 mmHg /sec

### **MEASUREMENT RANGE**

- Systolic: +60 to +260 mmHg
- Diastolic: +25 to +240 mmHg
- Pulse: 40 bpm to 199 bpm

### **MEASUREMENT METHOD**

- Oscillometry

### **OVER PRESSURE LIMIT**

- 305 +/-mmHg

### **EQUIPMENT**

- Class II medical device manufactured in FDA-regulated facility

### **ACCURACY**

- Compliant with standards



### **DIGITAL DISPLAY**

- Systolic & diastolic pressure
- Pulse rate

### **INPUT POWER REQUIREMENTS**

- 117VAC,60Hz
- 2.5 amps
- 3-wire grounded outlet



(800) 800.1787

AmericanMedicalScreening.com